

Analysis of Skeletal Remains Using Deep mtDNA Sequencing

MiSeq[®] system aids in identification of skeletal remains by sequencing mtDNA hypervariable regions.

Introduction

Forensic researchers face significant challenges in identifying unknown remains found in mass graves, such as those located in the war-torn countries of Vietnam, Kosovo, and Somalia. Dental records and circumstantial evidence are often inconclusive. DNA collected from skeletal samples is typically highly degraded, contaminated, and present in small quantities. The poor quality of these samples limits the success of conventional short tandem repeat (STR) typing techniques. In contrast, mitochondrial DNA (mtDNA) sequencing is quite useful in analyzing these samples, offering the following advantages:

- It is abundant and present in hundreds of copies per cell compared to nuclear DNA.
- It has a low mutation rate.
- It is maternally inherited and does not undergo sexual recombination, enabling maternal relationships to be established across generations.

For years, forensic laboratories have used capillary electrophoresis (CE)-based Sanger sequencing to perform mtDNA analysis. While the full mtDNA genome is approximately 16 kb, the hypervariable I and II (HVI and HVII) regions are usually the targets of mtDNA sequencing for human identification. These regions are located within a 1.1 kb region called the “D-loop”. Capable of providing deeper coverage of specific regions of interest, next-generation sequencing (NGS) systems enable mtDNA sequencing to be performed at a high throughput, using simpler workflows.

The Illumina MiSeq system is a fully integrated NGS instrument, providing a rapid workflow and streamlined sample preparation that is ideal for sequencing small genomes and amplicon samples, such as the D-loop region of mtDNA. This application note describes the use of the MiSeq system for deep sequencing of the HVI and HVII regions in mtDNA and how this technique is being used to aid the identification of human skeletal remains found in mass graves in Vietnam (Figure 1).

mtDNA Sequencing of Bone Fragments


Forensic scientists used the MiSeq system to sequence bone fragments uncovered at a mass grave site in Vietnam after DNA analysis of the samples using CE-based sequencing failed. The process began with the extraction of total DNA from blood and bone samples using established methods. Four individual amplicons spanning the HVI and HVII regions were amplified from mtDNA using PCR primers specific to the D-loop/Control region^{1,2} (Figure 2). Library preparation was performed using the Nextera[®] XT DNA Sample Preparation Kit. The samples were then sequenced on the MiSeq system using

Figure 1: Recovery of Human Remains



MiSeq technology is helping to identify human skeletal remains such as these, recovered in Vietnam. There are an estimated 300,000 unidentified human skeletal remains lying in mass graves across Vietnam.

Figure 2: NGS mtDNA Sequencing Workflow



Specific amplicons spanning the HVI and HVII regions were sequenced on the MiSeq system using a 2 x 150 bp run protocol.

a 2 × 150 bp run protocol. Raw reads were aligned to the mtDNA reference using BWA³, and variants were screened against the mtDNA revised Cambridge reference sequence (rCRS⁴) and called using the Genome Analysis Toolkit (GATK^{5,6}). The frequencies of each variant were calculated and compared to establish shared maternal lineage.

Analysis and Results


Concordant Sequencing Results of HVI/II Regions

Using this method, ultra-deep sequencing coverage (>10,000 reads) was reproducibly obtained for all samples. All variants previously reported in control DNA 9947A (Promega) were identified at high frequencies (Figure 3). There was concordance in variant positions and frequencies between technical replicates. Results for blood samples were concordant with those previously analyzed using CE-based Sanger sequencing.

Human Identification from Buried Remains


Total DNA was extracted from a bone sample (995A) recovered from a grave at the Southern Battlefield Cemetery in Vietnam. A corresponding blood sample (995B) was also drawn from the presumptive brother of the interred individual and sequenced using Illumina sequencing technology. The variant positions and frequencies in the HVI/II regions were compared (Figure 4). Results indicated that the two samples may share maternal lineage, supporting the hypothesis of kinship.

Figure 3: SNP Results Using MiSeq System Show Concordance



Raw frequencies are shown for each SNP across the D-loop. Blood and control results were concordant with samples previously sequenced with Sanger sequencing.

Figure 4: Illumina Sequencing of Bone Samples Confirms Putative Relationships



Note: Two pairs of amplicons cover each of the HV regions, with overlap of approximately 240 bp within HVI and 110 bp overlap within HVII.

DNA from a bone sample recovered from a grave (995A) and a blood sample (995B) were compared. The high degree of concordance in identified SNPs and profile characteristics suggests the two individuals are brothers.

